



## ICEROSS® SEAL-IN® X TF

Make your move

Iceross Seal-In TF liners are a trusted choice of above-knee amputees seeking a comfortable expulsion solution. However, some users have irregular limb shapes, or other anomalies, that conflict with the location of the seal.

Introducing Seal-In X TF, with proven Volume Adaptive Blade technology and a moveable seal ring that can be optimally positioned to maximize comfort for each unique individual. Seal-In X TF provides the personalized fit required to satisfy the varied needs of your above-knee amputees.



*Seal-In X TF features a moveable seal ring with Volume Adaptive Blades that change orientation to accommodate unique limb shapes and volume fluctuations.*

# ICEROSS® SEAL-IN® X TF

User Information	
Amputation Level:	Above-Knee (Transfemoral)
Impact Level:	All
Liner Information	
Sizes Standard:	25, 26, 28, 30, 32, 34, 36, 38, 40, 45, 50, 55
Sizes Conical:	23, 25, 26, 28, 30, 32, 34, 36, 38, 40
Size Seal Ring:	26, 28, 30, 32, 35, 38, 41, 44, 47, 51, 55, 60, 65, 70
Profile:	Conical & Standard



SEAL-IN



CONICAL



ACTIVE SKINCARE



FABRIC COVER



SILKEN INNER SURFACE

## ICEROSS SEAL-IN X TF LINER AND RING

A complete X TF system includes either a standard or conical profile liner and one seal. Replacement seals may be ordered separately.

Part#	Description
I-8532XX	Standard Profile Liner
I-8632XX	Conical Profile Liner
I-SXF0XX	Iceross Seal-In X TF Ring
I-SXF1XX	Iceross Seal-In X TF Replacement Ring

Ring size	Liner Size: Standard ■ Conical ●													
	23	25	26	28	30	32	34	36	38	40	45	50	55	
26	●	■												
28	●	■	■	■										
30	●	●	■	■	■									
32	●	●	●	■	■	■								
35	●	●	●	●	●	■	■							
38	●	●	●	●	●	●	■	■	■					
41	●	●	●	●	●	●	■	■	■	■				
44			●	●	●	●	●	■	■	■				
47				●	●	●	●	●	■	■	■			
51					●	●	●	●	●	●	■	■		
55						●	●	●	●	●		■	■	
60							●	●	●	●			■	
65								●	●	●				
70									●	●				

## MAXIMUM PERFORMANCE

New, longer length liner with a durable textile for stabilizing soft tissues and a redesigned distal end for a more intimate fit and enhanced proprioception.




## THINNER DISTAL END WITH FLEXIBLE COVER CLOSURE

Designed to better conform to limb shapes, Seal-In X TF helps eliminate air pockets when donning, thereby enhancing proprioception.

## MAXIMUM FLEXIBILITY

Choose a seal ring size that best suits the individual's physiology and move it to the perfect position for optimal suspension.



Measure for seal ring size (cm) directly on skin at lower edge of ring placement

FOLLOW ÖSSUR ON



TEL (800) 233-6263  
FAX (800) 831-3160  
WWW.OSSUR.COM

