

# Alpha<sup>®</sup> Hybrid Sleeve

## Product Highlights

- Provides extra protection at socket brim line
- Features Hybrid gel
- Flexible fabric


Featuring a highly flexible fabric for greater comfort, flexibility and increased range of motion, the Alpha Hybrid Sleeve may be used with locking suspension, suction suspension, or with LimbLogic<sup>®</sup> for creating an airtight seal. The Alpha Hybrid Sleeve uses Hybrid Gel which provides a balance of greater durability with comfort. For even more durability, the Alpha Hybrid Sleeve features a contoured polyurethane pad on the interior fabric cuff. This pad, which covers the typical socket brim, helps reduce damage to the seal due to impacts in the brim area.



✓ Transtibial

L-Code Suggestion: L5685

## Sleeve Sizing:



\*Hybrid Sleeve sizes 3-8 have staggered ranges of measurements.

## Part Numbers:

Sleeve Size	Buff Item #	Black Item #
Size 2	ALS-52B	ALS-52K
Size 3	ALS-53B	ALS-53K
Size 4	ALS-54B	ALS-54K
Size 5	ALS-55B	ALS-55K
Size 6	ALS-56B	ALS-56K
Size 7	ALS-57B	ALS-57K
Size 8	ALS-58B	ALS-58K

## Alpha Hybrid Sleeve: Specifications

Product Color: Buff, Black

Sleeve Thickness: 3 mm

Warranty: 3 months from date of patient fitting