

# Alpha<sup>®</sup> SmartTemp

Silicone AK Liners

## Product Highlights

- Liner utilizes Outlast<sup>®</sup> technology to absorb and store heat from a residual limb
- Delays the onset of sweat
- One-way select stretch fabric controls pistoning

The Alpha SmartTemp Silicone AK Liner featuring Outlast<sup>®</sup> combines silicone with the heat management technology of Outlast<sup>®</sup>. This revolutionary liner regulates heat by responding to fluctuations in skin temperature. By releasing stored heat as the body cools down, the liner stabilizes skin temperature to keep an amputee comfortable all day.


The Alpha SmartTemp Silicone AK Liner utilizes one-way stretch Select Fabric to control pistoning without the need of a distal matrix while allowing circumferential stretch for easy donning. The Select Fabric provides gentle tissue compression without placing uncomfortable pressure on the residual limb.



☑ Transfemoral

L-Code Suggestion: L5679 (cushion), L5973 (locking)

### Liner Sizing:



### Part Numbers:

Size/Profile	Locking Item #	Cushion Item #
Symmetrical AK		
Medium Plus 1	T442-1694-1	T440-1694-1
Medium Plus 2	T442-1694-2	T440-1694-2
Large 1	T442-1696-1	T440-1696-1
Large 2	T442-1696-2	T440-1696-2
Large Plus	T442-1698-1	T440-1698-1
Extra Large	T442-1697-1	T440-1697-1