

Alpha[®] Classic

MP Liners

Product Highlights

- Same price for locking or cushion
- Powerhouse core of the Alpha Family of liners
- Heat moldable


The most popular and best-selling Alpha Liner is the Alpha Classic MP. This liner features Alpha Classic Gel, 6 mm gel thickness, and the Original Fabric in green/gray.



✔ Transtibial

L-Code Suggestion: L5679 (cushion), L5973 (locking)

Liner Sizing:



Part Numbers:

Size	Cushion Item #	Locking Item #
Small	ALC-5063-E	ALL-5063-E
Medium	ALC-5060-E	ALL-5060-E
Medium Plus	ALC-5064-E	ALL-5064-E
Large	ALC-5066-E	ALL-5066-E
Large Plus	ALC-5068-E	ALL-5068-E
Extra Large	ALC-5067-E	ALL-5067-E