



# OVERLAY™



The Overlay is a breathable fabric sleeve with an integrated air expansion system that the wearer can adjust for maximum, personalized comfort.

## Device description

The Overlay™ system incorporates a set of air cells that enable custom dimensional adjustment for users. This system creates a uniform tightness in the socket, reducing pressure points and compensating for the residual limb volume variation when using a prosthesis. This product can be combined with most of the existing liners on the market.

## How to measure

Full product number example OV23-LG-L



C = Circumference at 4cm from the distal end with the liner on.

L = Length from the popliteal fossa to the tibial distal end with the liner on and the leg at 90°

O = Orientation of the pump to the user's preference.

CIRCUMFERENCE	
C (cm)	Circumference code
23,5 to 28,5	OV23
28,5 to 35,5	OV28
35,5 to 43,5	OV35

LENGTH	
L (cm)	Length code
14 to 19	SH (Short)
20+	LG (Long)

ORIENTATION	
O	Orientation
The pump is to the right when the product is worn	R
The pump is to the left when the product is worn	L

### Available sizes

- OV23-LG-L
- OV23-LG-R
- OV23-SH-L
- OV23-SH-R
- OV28-LG-L
- OV28-LG-R
- OV28-SH-L
- OV28-SH-R
- OV35-LG-L
- OV35-LG-R
- OV35-SH-L
- OV35-SH-R

